

Team Quadcopter: Code

Arduino Code

```
// Description: This is an Arduino program to measure the value of a resistor
in ohms.

// Include bluetooth libraries
#include <SPI.h>
#include <boards.h>
#include <RBL_nRF8001.h>
#include <services.h>

// Initialize global variables
int vinSupplier= A0; // The 5V input voltage pin
int vinReader = A2; // The pin that reads the exact input voltage
int voutReader = A1; // The pin that reads the voltage accross the resistor
being measured
int triggerPin = A3; // The pin that listens to whether the button used to
initiate the measurement is pressed
int ledIndicator = A4; // Pin used to drive LED to alert the user a measurement
is in progress
int voutRaw= 0; // A variable to hold the ADC measurement of the voltage
accross measured resistor
int vinRaw = 0; // A variable to hold the ADC measurement of the input voltage
float Vin= 0; // vinRaw mapped to an actual voltage value
float Vout= 0; // voutRaw mapped to an actual voltage value
float Rknown= 14840;//14840;//991; // the value of the known resistor
float Resistance= 0; // the calculation of resistance
char charRes[10]; // The character array representation of the resistance for
ble transmission
int triggerPinVoltage = 0; // the ADC measurement of the switch voltage

boolean inprogress = false; // a boolean to let us know whether a test is
currently in progress
boolean canceled = false;
int avgLength = 1000; // The number of measurements to be averaged
float average(float a[]); // An averaging function

void setup()
{
    // Set the vinSupplier pin to be an output
    pinMode(vinSupplier, OUTPUT);
    pinMode(ledIndicator, OUTPUT);
}
```

```

    // Begin serial communication
    Serial.begin(9600);

    // Init. and start BLE library.
    ble_begin();
}

void loop()
{
    // Check to see if the push button is pressed to initiate a measurement
    triggerPinVoltage = analogRead(triggerPin);

    // Only measure if we have received a ble command to do so
    if( ble_available() ){
        analogWrite(ledIndicator, 255);
        Serial.println("Ble Available");

        // Exhaust the ble data buffer so we have nothing left in there to cause
        // a second test to run without user prompt
        while ( ble_available() )
            ble_read();

        // Let's take a whole bunch of measurements and average
        // them to increase precision
        //float resistanceArray[avgLength];
        float sum = 0;
        int i;
        for(i=0; i<avgLength && canceled == false; i++){

            if(ble_available()){
                canceled = true;
                Serial.println("test canceled");
                while ( ble_available() )
                    ble_read();
            }
            ble_do_events();
            // Set inprogress to true to indicate a test has begun due to
            // a button press
            inprogress = true;

            // Drive the input to the ohm meter with a 5V input
            analogWrite(vinSupplier, 255);

            // Read the voltage across the measured resistor
            voutRaw= analogRead(voutReader);

            // Read the voltage being supplied by vinSupplier
            vinRaw = analogRead(vinReader);

```

```

    // Turn off the input voltage
    analogWrite(vinSupplier, 0);

    // Map the ADC measured voltages to actual voltage values
    Vin = (vinRaw/1024.0)*5.0;
    Vout = (voutRaw/1024.0)*5.0;

    // Compute the current flowing through the known resistor
    float current = (Vin - Vout)/Rknown;

    // Determine the resistance of the resistor we are measuring
    Resistance = Vout/current;
    sum = sum + Resistance;
    Serial.println(Resistance);
    //Serial.println(',');

    // Turn off the input voltage
    //analogWrite(vinSupplier, 0);

    //Add the resistance value to our averaging array
    //resistanceArray[i] = Resistance;
    //delay(1000);

} // End for loop

canceled = false;

//float finalResistance = average(resistanceArray);
float finalResistance2 = sum/avgLength;
Serial.println("Final Resistance");
// Serial.println(finalResistance);
Serial.println(finalResistance2);

// Convert float resistance to char array so that we
// can transmit over ble
dtostrf(finalResistance2, 4, 3, charRes);

for(int i=0;i<sizeof(charRes);i++)
{
    // Write to ble and to console
    ble_write(charRes[i]);
    Serial.print(charRes[i]);
}

ble_write('\n'); //the new line character let's the phone know the

```

```

message is over
    analogWrite(ledIndicator, 0);

}
else if(triggerPinVoltage > 10 && inprogress){
    // Do not do anything if the button has remained depressed since
    // the latest resistance measurement
    //Serial.println("Idling");

}
else{
    // This means the button has be unpresed so we can finally
    // say this resistance measurement is over and can start a new one
    // once the button is pressed again.
    inprogress = false;
}

// Communicate changes
ble_do_events();
//delay(1000);

}

float average(float a[]){
    int i;
    float avg, sum=0.0;
    for(i=0;i<avgLength;++i){
        sum+=a[i];
    }
    avg =(sum/avgLength);
    return avg;
}

```

App Code

```

angular.module('starter', ['ionic'])

.run(function($ionicPlatform, $templateCache, $rootScope, $interval, $timeout)
{

    // Load a random value to populate the knob
    $rootScope.currentValue = 41;

```

```

$scopeScope.weight = 0;

// Create a global ble connected status variable for use with ng-show/ng-hide
$scopeScope.connected = false;

$ionicPlatform.ready(function() {

    // Hide the accessory bar by default (remove this to show the accessory bar
above the keyboard
    // for form inputs)
    if(window.cordova && window.cordova.plugins.Keyboard) {
        cordova.plugins.Keyboard.hideKeyboardAccessoryBar(true);
    }
    if(window.StatusBar) {
        StatusBar.styleDefault();
    }

    // monitor ble connectivity status every second
    $scopeScope.monitorConnectivity = function(){
        $interval(function(){
            bluetoothSerial.isConnected(
                function() {
                    console.log("bluetooth connected");
                    $scopeScope.connected = true;
                },
                function() {
                    $scopeScope.connected = false;

                    // If we are not already trying to connect, attempt to
reconnect
                    if(!$scopeScope.connecting){

                        $timeout(function(){
                            console.log("attempting to RECONNECT")
                            $scopeScope.connecting = true;
                            $scopeScope.initiateConnection();

                        }, 1000)

                    }
                    console.log("Bluetooth is *not* connected");
                }
            );
        }, 1000)
    }();

    $scopeScope.initiateConnection = function(){
        console.log("running initiateConnection..")
    }

```

```

    // Let's scan the environment for a ble device with the name "BLE Shield"
    bluetoothSerial.list(function(devices) {
        console.log("in success callback")

        if(devices.length === 0)
            return $rootScope.connecting = false;

        var bleShield = _.findWhere(devices, {name: "BLE Shield"});

        // JSON output of bleShield looks like this:
        // [{"id":"55174456-779D-D60E-82D4-EA927560790C","name":"BLE
Shield","uuid":"55174456-779D-D60E-82D4-EA927560790C"}]

        // Now that we have found the bleShield and have its id, lets
connect!
        bluetoothSerial.connect(bleShield.uuid, $rootScope.connnectSuccess,
function(){ $rootScope.connecting = false; });
        }, function(){
            console.log("in fail callback")
        });
    }

    $rootScope.gotMessage = function(data){
        // This function gets called when we receive the resistance measurement
from the arduino
        console.log("got gotMessage called")
        console.log(data)
        $rootScope.newResistance = parseFloat(data);
    }

    $rootScope.connnectSuccess = function (){
        console.log("We connected via bluetooth");
        console.log("attempting to subscribez")
        $rootScope.connecting = false;
        // $rootScope.connected = true;
        // $scope.connected = true;
        console.log("connected: " + $rootScope.connected)

        bluetoothSerial.subscribe("\n", $rootScope.gotMessage, function(){
            console.log("got an error..");
        });
    }

    });
})

```

```

.config(function($stateProvider, $urlRouterProvider) {
  $urlRouterProvider.otherwise('/')

  $stateProvider
    .state('home', {
      url: '/',
      //template: '<p style="color: white;">Hello, world!</p>'
      templateUrl: 'home.html',
      controller: 'homeCtrl'
    })

    .state('enterWeight', {
      url: '/enterWeight',
      templateUrl: 'enter-weight.html',
      controller: 'enterWeightCtrl'
    })

    .state('measuring', {
      url: '/measuring',
      templateUrl: 'measuring.html',
      controller: 'measuringCtrl'
    })
  })

.controller('homeCtrl', function($scope, $rootScope, $timeout) {

  $timeout(function(){

    if(!$rootScope.connected)
      $rootScope.initiateConnection();

  }, 1000)

  function connectFailure(){
    console.log("Could not connect");
  }

  function gotMessageError(error){
    // If there is an error receiving the resistance measurement from the
    arduino, this function
    // gets called
    console.log("error")
  }

  $timeout(function(){

    // Run the code to initiate our body fat percentage knob (this should

```

```

ideally be within a directive)
    initiateKnob();
    initiateLoader($rootScope.currentValue);

    // Set knobInitiated to true so we can unhide the now not-so-ugly knob
    $scope.knobInitiated = true;
  })

  })

.controller('enterWeightCtrl', function($scope, $rootScope, $timeout,
$interval, $state) {
    $scope.fullWeight = 0;

    $scope.storeWeight = function(){
        $rootScope.weight = parseFloat($scope.fullWeight);
        $state.go("measuring");
    }
    //console.log($rootScope.currentValue)
  })

.controller('measuringCtrl', function($scope, $rootScope, $timeout, $interval,
$state) {
    console.log($rootScope.currentValue)

    delete $rootScope.newResistance;

    bluetoothSerial.write("start_test", function(){
        console.log("started_test");
    }, function(){
        console.log("failed to start_test");
    });

    $scope.progressValue = 0;

    incrementProgress = $interval(function(){
        console.log("interval called")
        if($scope.progressValue < 99 && typeof $rootScope.newResistance ===
"undefined")
            $scope.progressValue++;

        else if(typeof $rootScope.newResistance === "undefined"){
            // wait to get resistance
        }
        else {
            $scope.progressValue = 100;
            // cancel the interval since it appears to be global accross all

```



```

controllers
    $interval.cancel(incrementProgress);

    console.log("about to compute regression")
    // Define our leanMass regression coefficients
    // var B0 = -13.36117;
    // var B1 = 0.0085998;
    // var B2 = 0.5964701;

    var B0 = 1.207298;
    var B1 = -0.0047503;
    var B2 = 0.5542703 ;

    // Define our waterMass regression coefficients
    var C0 = 14.35698;
    var C1 = -0.0079816;
    var C2 = 0.3226745;

    // $rootScope.newResistance = 600;
    console.log("$rootScope.weight: " + $rootScope.weight);
    // $rootScope.weight = 130;
    // Let's compute the regression jquery.classyloader.js
    var leanMass = B0 + B1*$rootScope.newResistance + B2*$rootScope.weight;
    console.log("leanMass: " + leanMass)
    var waterMass = 0.35*$rootScope.weight; // C0 + C1*$rootScope.newResistance
+ C2*$rootScope.weight;
    console.log("waterMass: " + waterMass)
    var fatMass = $rootScope.weight - leanMass - waterMass;
    console.log("fatMass: " + fatMass)

    // set a random value for the new body fat percentage from (0-100)
    $rootScope.currentValue = Math.round(100*(fatMass/($rootScope.weight -
waterMass))); //Math.floor(Math.random() * 100) + 1;
    // $rootScope.currentValue = 21.9;
    console.log("fatPercentage: " + $rootScope.currentValue);
    delete $rootScope.newResistance;
    // Go back to the home state to display the new bodyfat percentage
    $state.go('home');
}
}, 80)

$scope.cancelMeasurement = function(){
    $interval.cancel(incrementProgress);
    console.log("cancelMeasurement called")
    bluetoothSerial.write("cancel_test", function(){
        console.log("canceled test");
        $state.go('home');
    }, function(){

```

```

        console.log("failed to cancel_test");
    });

}

})

function initiateKnob(){

    // This function is used to initialize the jQuery knob used to display the
    body fat percentage
    // to the user in the home state

    $(function($) {
        $(".knob").knob({
            change : function (value) {
                //console.log("change : " + value);
            },
            release : function (value) {
                //console.log(this.$.attr('value'));
                console.log("release : " + value);
            },
            cancel : function () {
                console.log("cancel : ", this);
            },
            /*format : function (value) {
                return value + '%';
            },*/
            draw : function () {
                // "tron" case
                if(this.$.data('skin') == 'tron') {
                    this.cursorExt = 0.3;
                    var a = this.arc(this.cv) // Arc
                        , pa // Previous arc
                        , r = 1;
                    this.g.lineWidth = this.lineWidth;
                    if (this.o.displayPrevious) {
                        pa = this.arc(this.v);
                        this.g.beginPath();
                        this.g.strokeStyle = this.pColor;
                        this.g.arc(this.xy, this.xy, this.radius -
this.lineWidth, pa.s, pa.e, pa.d);
                        this.g.stroke();
                    }
                    this.g.beginPath();
                    this.g.strokeStyle = r ? this.o.fgColor :
this.fgColor ;
                    this.g.arc(this.xy, this.xy, this.radius -

```

```

this.lineWidth, a.s, a.e, a.d);
        this.g.stroke();
        this.g.lineWidth = 2;
        this.g.beginPath();
        this.g.strokeStyle = this.o.fgColor;
        this.g.arc( this.xy, this.xy, this.radius -
this.lineWidth + 1 + this.lineWidth * 2 / 3, 0, 2 * Math.PI, false);
        this.g.stroke();
        return false;
    }
}
});
// Example of infinite knob, iPod click wheel
var v, up=0,down=0,i=0
    , $idir = $("div.idir")
    , $ival = $("div.ival")
    , incr = function() { i++; $idir.show().html("+").fadeOut();
$ival.html(i); }
    , decr = function() { i--; $idir.show().html("-").fadeOut();
$ival.html(i); };
    $("input.infinite").knob(
        {
            min : 0
            , max : 20
            , stopper : false
            , change : function () {
                if(v > this.cv){
                    if(up){
                        decr();
                        up=0;
                    }else{up=1;down=0;}
                } else {
                    if(v < this.cv){
                        if(down){
                            incr();
                            down=0;
                        }else{down=1;up=0;}
                    }
                }
                v = this.cv;
            }
        }
    );
});

}

function initiateLoader(pct){

```

```
$('.loader').ClassyLoader({
  speed: 50,
  diameter: 80,
  fontSize: '30px',
  fontFamily: 'Arial',
  fontColor: 'rgb(53, 188, 228)',
  lineColor: 'rgb(53, 188, 228)',
  remainingLineColor: 'rgba(73, 125, 164, 0.1)',
  percentage: pct,
  lineWidth: 20,
  start: 'top',
  //remainingLineColor: 'rgba(200,200,200,0.1)'
});
}
```